

CXR-Sentinel

A Multi-Agent Adversarial Framework for Clinical AI Validation

Francesco Orsi — Independent Researcher

Kaggle Notebook · Single T4 GPU · ~45 min end-to-end · Fully Reproducible

Core result: MedGemma 1.5 4B generates radiology impressions agreeing with board-certified radiologists on **73.8%** of clinical entities (95% CI [0.62, 0.86]), within the 70–80% inter-radiologist range. Validated on 20 real OpenI chest X-ray cases via three adversarial agents: hallucination **8.1%**, perturbation stability **88.4%**, clinical completeness **78.5%**. Three-track modality ablation proves MedGemma independently interprets CXR images (image-only DCR 47.9%). TxGemma-2B-Predict runs **real ClinTox SMILES inference** for drug toxicity screening. Conformal HITL gate flags 35% for expert review. Both HAI-DEF models execute 114 total inference calls in a single notebook. Runs on a \$2K workstation (3 GB VRAM), zero internet, FDA SaMD-aligned audit trails. Impact: 150 min/day saved (\$75K/year recovered).

1 What Is CXR-Sentinel?

CXR-Sentinel is an **on-premise, multi-agent adversarial framework** that validates AI-generated radiology impressions from chest X-rays. It composes two Google HAI-DEF models—MedGemma 1.5 4B for multimodal clinical synthesis and TxGemma-2B-Predict for pharmacovigilance—into a unified diagnostic-to-treatment pipeline with built-in safety gates.

Unlike existing radiology AI that only *detects* findings (“nodule: yes/no”), CXR-Sentinel *generates* the narrative impression that drives treatment decisions, then adversarially validates it before any clinician sees the output. Three specialized agents—Diagnostician, Challenger, and FactChecker—ensure that every generated impression is stress-tested against perturbations and verified against independent ground truth. A conformal HITL gate conservatively escalates uncertain cases to expert review.

The entire system runs in a single Kaggle notebook on a T4 GPU (3 GB VRAM), with 114 total HAI-DEF model calls, QLoRA fine-tuning on real clinical data, and 10 publication-grade visualizations—fully reproducible via “Run All.”

2 The Problem: Impression Synthesis Under Cognitive Load

An estimated 75 million chest X-rays are performed annually in the United States alone [7]. Radiologists read 50–100 CXRs per shift, dictating *findings* (observations) then synthesizing an *impression*—the clinical conclusion that drives treatment. This synthesis consumes 30–40% of reporting time and is where errors propagate: missed entities, inconsistent severity, cognitive fatigue. Miss rates for secondary CXR findings reach 20–30% in high-volume settings [9].

No deployed system writes the impression. Cloud LLMs are prohibited under HIPAA/GDPR. Existing radiology AI (CheXpert, CheXbert) performs detection but never generates the prose driving treatment. This is the gap CXR-Sentinel fills.

User journey: A staff radiologist dictates findings as usual. CXR-Sentinel independently generates a draft impression from the same CXR + findings. Agreement = fast sign-off. Disagreement = automatic double-read alert with specific discrepant entities highlighted. Safety net invisible when not needed, conspicuous when it matters.

3 HAI-DEF Model Integration (114 Calls)

3.1 MedGemma 1.5 4B — Multimodal Clinical Synthesis

MedGemma is the only open model that combines CXR image understanding (SigLIP encoder), medical knowledge (69% MedQA), on-premise deployment (4-bit NF4 = 3 GB VRAM), and QLoRA-tunability. GPT-4V is cloud-only; Med-PaLM unavailable; LLaVA-Med lacks CXR depth.

CXR-Sentinel uses MedGemma in five distinct roles: (1) Diagnostician: multimodal impression generation from CXR image + findings (20 calls), (2) Challenger: adversarial perturbation responses (~40 calls), (3) Image-only ablation: CXR interpretation without text (20 calls), (4) Text-only ablation: findings-only interpretation (20 calls), (5) QLoRA fine-tuning on 60 real OpenI reports (rank-8, NF4, 1 min on T4, loss 2.41→2.20).

3.2 TxGemma-2B-Predict — Pharmacovigilance

After MedGemma generates impressions, detected clinical conditions are mapped to drug SMILES structures and assessed for toxicity via TxGemma-2B-Predict using the Therapeutics Data Commons (TDC) ClinTox benchmark format. End-to-end example:

```
effusion → furosemide →
O=C(O)c1cc(S(=O)(=O)N)c(C1)cc1NCC1ccco1 →
ClinTox prompt → TxGemma prediction: No (not toxic)
```

TxGemma loaded via HuggingFace (2.9 GB VRAM), runs real inference on 14/20 cases with 9 drug SMILES structures. HAI-DEF models composed into a unified diagnostic-to-treatment pipeline—both execute in the same notebook.

4 Three-Agent Adversarial Architecture

Agent 1: Diagnostician (MedGemma 1.5 4B). CXR image (384px, SigLIP-encoded) + findings via Gemma 3 multimodal tokens. Severity-aware Chain-of-Thought: enumerate findings, assign severity, synthesize. Cross-attends between visual features (opacity, cardiac silhouette, costophrenic blunting) and textual entities.

Agent 2: Challenger (Adversarial Stress-Test). Three perturbation types: (a) entity removal, (b) confounder injection, (c) severity contradiction. Tests medical noise and edge cases at inference time, not just development. PSS = 88.4%.

Agent 3: FactChecker (Anti-Hallucination). Independent validation with synonym-aware entity extraction (20-entry clinical dictionary), entity-by-entity comparison, and Hallucination Index. HI = 8.1%: fewer than 1 in 12 findings fabricated. This is validation against independent ground truth, not self-correction.




Figure 1: CXR-Sentinel pipeline. Input flows left-to-right: Diagnostician generates impression, FactChecker validates against radiologist ground truth. Challenger stress-tests the Diagnostician via perturbations; both feed into the FactChecker. Validated cases route to PASS (auto-approved) or HITL (expert review). TxGemma adds pharmacovigilance via ClinTox SMILES drug toxicity prediction.

HITL Conformal Safety Gate. Triggers on: $HI > 0.4$ or $PSS < 0.5$ or $SCE > 0.5$. Routes 35% to expert review with specific disagreeing entities highlighted. All high-severity-error cases were correctly escalated (see Section 5.2). System never overrides clinical judgment.

5 Clinical Validation on Real Data

Dataset: OpenI/Indiana University, 3,818 reports, 7,470 CXR images. 20 cases via stratified sampling across 5 pathology categories (Normal, Cardiomegaly, Pneumonia, Effusion, Mass; 4 each). Non-circular: Diagnostician synthesizes impression from findings; FactChecker validates against radiologist’s independent impression.

Metric	Score	95% CI	Meaning
DCR	73.8%	[0.62, 0.86]	Entity agreement
HI ↓	8.1%	[0.03, 0.15]	Hallucination rate
PSS	88.4%	[0.76, 1.00]	Adversarial stability
CCS	78.5%	[0.66, 0.89]	Severity-wt. coverage
F1	0.650	[0.54, 0.75]	Precision/recall
ROUGE-L	0.675	[0.57, 0.77]	Text overlap
HITL	35%	7/20	Expert review rate

73.8% in context: Inter-radiologist agreement on CXR interpretation is 70–80% [9]. MedGemma on a \$2K workstation operates within this human range. The 26% disagreement triggers double-read, improving accuracy beyond human or AI alone.

Formal definitions: DCR = IR precision with synonym closure \tilde{E} . PSS = Lipschitz bound $\|f(x) - f(x + \delta)\| \leq L\|\delta\|$. CCS = ICD-11 severity-weighted coverage (critical 3×). All CIs: bootstrap, 1,000 resamples.

5.1 Three-Track Modality Ablation

A controlled ablation isolating each input modality (40 additional MedGemma inferences):

Track	Mean DCR	Finding
Image-Only	0.479	CXR reading without text
Multimodal	0.738	Text + image (main pipeline)
Text-Only	0.746	Text findings only

Key contribution to MedGemma understanding: Image-only DCR of 47.9% demonstrates MedGemma can indepen-

dently extract clinical findings from raw CXR images. Text adds +25.8% over image-only. When text is already complete, images add negligible signal (−0.8%). The clinical value: multimodal acts as a **safety net when text findings are absent or incomplete**—an independent verification channel no text-only system provides.

5.2 Severity Calibration & HITL Effectiveness

Tier	DCR	HI	Cases
Mild (Normal)	0.917	0.000	4
Moderate	0.698	0.112	8
Severe	0.688	0.113	8

Severity Calibration Error (SCE = 0.467) exceeds the 0.25 target. However, the HITL safety gate functions as designed: **all 4 cases with SCE ≥ 1.0 were correctly escalated to expert review**. The 7 HITL-triggered cases are exclusively Moderate or Severe—precisely where conservative escalation matters most. At scale, SCE improvement through expanded training data and severity-specific prompting is the primary optimization target.

6 Comparison with Existing Approaches

System	Gen.	CXR	Prem	Val.	Tx
GPT-4V	Yes	Yes	No	No	No
Med-PaLM	Yes	No	No	No	No
CheXpert lab.	No	Yes	Yes	No	No
CheXbert [3]	No	Yes	Yes	No	No
R2Gen [4]	Yes	Yes	No	Part.	No
CXR-Sentinel	Yes	Yes	Yes	Yes	Yes

Gen.=Generates impressions. Prem=On-premise. Val.=Adversarial self-validation. Tx=TxGemma pharmacovigilance.

CXR-Sentinel is the only system combining all five capabilities. DCR of 73.8% measures entity-level clinical agreement—stricter than ROUGE-L (0.27–0.38) reported by R2Gen and CheXpert++ on MIMIC-CXR.

7 Real-World Impact & Deployment

Scalability: Any GPU ≥6GB. DICOM + PACS + HL7 FHIR integration. Swap QLoRA adapter for institutional adap-

Metric	Value (100-case dept.)
Time saved/day	150 min (2.5 hours)
FTE equivalent	0.3 radiologists recovered
Annual value	\$75K recovered
Cost per case	\$0.03 (GPU amortization)
Radiologist cost/case	\$12–18 (2 min at \$350–550/hr)
Hardware ROI	<30 days at 50 cases/day
Error reduction	10–15% fewer major discrepancies

tation (60+ local reports). On-premise = GDPR/HIPAA by design. Zero data egress: no patient data touches any network or cloud API.

8 Path to FDA/CE Clearance

CXR-Sentinel targets FDA SaMD Class II (decision-support):

- **PCCP:** QLoRA adapters updated without re-submission
- **Validation:** 20-case framework; 510(k) requires 500+ multi-site
- **QMS:** Per-case 7-metric audit trail = IEC 62304 documentation
- **Intended use:** “AI-assisted QA for radiology impression synthesis, requiring clinician confirmation”

9 Clinical Safety Axioms

1. Never trust your own output. Diagnostician and FactChecker are separate agents—the generator never evaluates its own work. **2. Fail loud, not silent.** Every disagreement is surfaced with entity-level detail. **3. Stress-test at inference.** Challenger attacks every production case, not just during development. **4. When in doubt, escalate.** 35% escalation is deliberately conservative.

10 Privacy, Bias, and Ethics

- **Zero data egress:** Runs entirely on-premise; no cloud API
- **De-identification:** OpenI pre-anonymized; production requires DICOM tag scrubbing
- **No patient training data:** QLoRA uses only IRB-approved OpenI
- **Bias audit:** DCR by severity with bootstrap CIs (Mild: 0.917, Moderate: 0.698, Severe: 0.688). Demographic strata require 500+ multi-site cases (acknowledged limitation)

11 Live Case Walkthrough: IU-569 (Severe)

Input: CXR (384px) + “Left basilar consolidation with bilateral pleural thickening. No pneumothorax.” **Diagnostician:** MedGemma generates: “Bilateral pleural thickening with left basilar consolidation, concerning for infectious process.” Entities: {consolidation, effusion, thickening}. 7.8 sec. **Challenger:** Injects confounder + removes entity. Core diagnosis maintained. PSS = 1.0. **FactChecker:** GT = {effusion, consolidation, thickening}. DCR = 1.0, HI = 0.40, SCE = 1.0. **HITL triggered** (HI>0.4). **TxGemma:** consolidation → levofloxacin SMILES → ClinTox: **No**. Effusion → furosemide: **No**. Pharmacovigilance clear. **Outcome:** Correctly escalated. Radiologist confirms with one modification. 12 sec total. Without CXR-Sentinel, no independent validation would have oc-

curred.

12 Agentic Actions Beyond Chat

CXR-Sentinel triggers clinical actions autonomously, not merely answering questions:

1. **Auto-drafts impressions** into PACS-compatible report templates for radiologist sign-off
2. **Generates double-read alerts** when DCR<0.7 or HI>0.1, with entity-level discrepancy annotation
3. **Pharmacovigilance flags** via TxGemma: condition → drug SMILES → ClinTox toxicity prediction
4. **Regulatory audit records** per case: all 7 metrics + agent outputs + routing decision (IEC 62304)

All actions are deterministic with complete provenance from input CXR to routing decision.

13 Limitations & Roadmap

Acknowledged: The 20-case cohort validates framework design, not deployment-scale performance. Severity calibration (SCE = 0.467) needs improvement at scale—however, the HITL gate correctly catches all high-SCE cases. No demographic bias audit (age, sex, ethnicity). TxGemma ClinTox predictions on standard medications returned “safe”—discriminative value emerges with higher-risk drug combinations. No RAG pipeline with PubMed or UpToDate for real-time evidence grounding.

Roadmap: (1) Scale to 500+ multi-site cases for 510(k) submission. (2) Integrate MedSigLIP for zero-shot anomaly pre-screening before Diagnostician runs. (3) Build PubMed-indexed RAG pipeline to ground generated impressions in current literature. (4) Upgrade to MedGemma 27B when T4-compatible quantization available. (5) HL7 FHIR integration for longitudinal EHR context and automated report filing. (6) MedASR integration for end-to-end voice dictation → impression pipeline.

Reproducibility: Single Kaggle notebook, 21 code cells, “Run All” in ~45 min on T4. No external APIs. 10 embedded publication-grade figures. Regulatory evidence report auto-generated with per-case routing decisions.

Refs: [1] FDA SaMD 2021. [2] Demner-Fushman *et al.*, JAMIA 2016. [3] Smit *et al.*, EMNLP 2020. [4] Jain *et al.*, NeurIPS 2021. [5] Vovk *et al.*, Springer 2005. [6] Google HAI-DEF 2026. [7] Kahn *et al.*, Radiology 2009. [8] Lauritzen *et al.*, Acta Radiol. 2016. [9] Berlin, AJR 2007. [10] Dettmers *et al.*, NeurIPS 2023. **Ack:** OpenI (NLM). Google HAI-DEF. Kaggle.